



▶ Semi-Flexibles Nephroskop

Semi-Flexible Nephroscope



grenzenlos im Detail! ▶▶




Semi-Flexibles Nephroskop

Semi-Flexible Nephroscope









Inkl. Reinigungsbürste für Arbeitskanal und Spülkanal
Incl. cleaning brush for forcep channel and jet channel

Technische Daten Technical Data

Standard Form / Standard Style Ø 7,5-9		Model: URS-NS
instrument size (Charr.):	Länge / length	220mm
	Außenrohr / outer tube dim.	7,5-9 12 Fr.
	Arbeitskanal / working channel	6,8 Charr.
	Blickrichtung / direction of view	5°
	Öffnungswinkel / field of view	85°
	Bildpunkte / pixel	40.000
	Artikel Nr. / item no.	240.00133.00

Zubehör / Accessories

Bild / Picture	Artikelnummer / Item No.	Beschreibung / Description	Maße / Dimensions	VE / PU
	210.00083.00	Set Ventilkappe + Hilfsdichtung dome valve + backup seal	ML 151 VMG + ML 151 VMQ	1 VE = 10 St. 1 PU = 10 pcs
	220.00025.00	Reinigungsbürste für Arbeitskanal cleaning brush for working channel	Bürste/brush Ø3,0/L=30mm, L=590mm	1St./pc
	220.00026.00	Reinigungsbürste für Spülkanal cleaning brush for flushing channel	Bürste/brush Ø3,0/L=10mm, L=590mm	1St./pc
	210.00088.00	Hahnreiber rechts tap cock - right		1St./pc
	210.00089.00	Hahnreiber links tap cock - left		1St./pc
	210.00090.00	Federkappe protection cap		1St./pc