

Hysteroskope

Hysteroscopes



Hysteroskope

Hysteroscopes



Technische Daten

Technical Data

Hysteroskop Ø 2.8 mm / Hysteroscope Ø 2.8 mm

Artikel Nr. / Item no.	HS28A.01	HS28B.01
Durchmesser / Diameter	Ø 2.8 mm	Ø 2.8 mm
Kompatibel / Compatible	Storz®	Storz®
Nutzlänge / Working length	302,5 mm	302,5 mm
Blickrichtung / Direction of view	0°	30°

Artikel Nr. / Item no.	HS28C.01	HS28F.01
Durchmesser / Diameter	Ø 2.8 mm	Ø 2.8 mm
Kompatibel / Compatible	Storz®	Storz®
Nutzlänge / Working length	302,5 mm	302,5 mm
Blickrichtung / Direction of view	70°	12°

Hysteroskop Ø 4.0 mm / Hysteroscope Ø 4.0 mm

Artikel Nr. / Item no.	HS40A.01	HS40B.01
Durchmesser / Diameter	Ø 4.0 mm	Ø 4.0 mm
Kompatibel / Compatible	Storz®	Storz®
Nutzlänge / Working length	302,5 mm	302,5 mm
Blickrichtung / Direction of view	0°	30°

Artikel Nr. / Item no.	HS40C.01	HS40F.01
Durchmesser / Diameter	Ø 4.0 mm	Ø 4.0 mm
Kompatibel / Compatible	Storz®	Storz®
Nutzlänge / Working length	302,5 mm	302,5 mm
Blickrichtung / Direction of view	70°	12°